


## Checking Your Neurostimulator Battery

Please check your battery on a regular basis.

When the screen shows **2.80V**, check your battery once a month.

When it shows **2.70V**, check your battery once a week.

If you see **ERI** on your remote,



Call the DBS clinic as soon as possible to set up surgery for battery replacement.

For more copies, go online at <http://vch.eduhealth.ca> or email [pthem@vch.ca](mailto:pthem@vch.ca) and quote Catalogue No. **FM.130.B17**

© Vancouver Coastal Health, October 2016

The information in this document is intended solely for the person to whom it was given by the health care team.

[www.vch.ca](http://www.vch.ca)